

ELTRA HELIOS CARBON SUFLUR DETERMINATOR

CONSUMABLES LIST

Alpha Part No.	Description
	Set of consumables for 1000 analyses consisting of: 1000 Ceramic crucibles (1 Set = 1000 pcs) 908g anhydrone 50g glass wool 500g sodium hydroxide
	or:
	Set of consumables for 1000 analyses consisting of: 1000 Ceramic crucibles (2 x 500 pcs = 1 Set) 908g anhydrone 50g glass wool 500g sodium hydroxide
AR-3818	Ceramic crucibles (2 x 500 pcs = 1 Set)
AR-171	Anhydrone
AR-177	Quartz wool
AR-265	Glass wool
AR-2174	Sodium hydroxide
AR-90300	Iron phosphate
AR-636	Vanadiumpentoxide
AR-90280	Tin accelerator
	Coal standard (select S value)*
AR-90800	Graphite (100%C)
AR-1034	Calcium carbonate (12%C)
AR-241	High vacuum grease

For a competitive quotation on your Consumables needs for all Eltra and Leco instruments please contact :

ALPHA CANADA INC.

Tel : 905.824.1198

Fax : 905.824.2199

Email : tanyar@alphacanadainc.com