

## ELTRA CS-530 CARBON SULFUR DETERMINATOR

### CONSUMABLES LIST

Alpha Part No.	Description
	Set of consumables for 1000 analyses consisting of: 1000 reusable boats 58x22x14mm 908g anhydrone 50g glass wool 500g sodium hydroxide
	or:
	Set of consumables for 1000 analyses consisting of: 1000 disposable boats 88x13x9mm 908g anhydrone 50g glass wool 500g sodium hydroxide
AR-8203	Ceramic boats, 58x22x14 mm
	Porcelain boats, 88x13x9 mm
AR-171	Anhydrone
AR-177	Quartz wool
AR-265	Glass wool
AR-2174	Sodium hydroxide
AR-90300	Iron phosphate
AR-636	Vanadium pentoxide
AR90280	Tin accelerator
AR-608	Halogen trap
	Coal standard (select S value)*
	Graphite (100%C)
AR-1034	Calcium carbonate (12%C)
AR-241	High vacuum grease

**For a competitive quotation on your Consumables needs for all Eltra and Leco instruments please contact :**

**ALPHA CANADA INC.**

**Tel : 905.824.1198**

**Fax : 905.824.2199**

**Email : [tanyar@alphacanadainc.com](mailto:tanyar@alphacanadainc.com)**